

Optimal Mattress Cover

Durable mattress topper with good fit & flexibility

Optimal Mattress Cover is a simple, durable and functional mattress topper with a good fit and flexibility. Optimal Mattress Cover is a removable and liquidproof hygiene cover for easy cleaning.

The liquidproof hygiene cover is manufactured in a four-way stretch⁽¹⁾⁽²⁾ fabric to reduce the risk of shear forces, and is vapour-permeable⁽³⁾ to lower the risk of skin maceration.



Teknisk specifikation	
Material	Polyurethane/polyester
Sizes	80×205 cm 85×200 cm 90×200 cm 105×200 cm 120×200 cm
Cleaning instructions	Wipe clean with detergent and/or disinfectant. Machine wash max. 95°C, tumble dry.
CE-marking	Registered and marked in accordance with MDR (EU) 2017/745

References

(1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.

(3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instructions for use prior to use.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2024. All rights reserved. Last update: 2024-01-22.